

# Rebound Tonometer. Model.SW-500

AN ECONOMICAL, ACCURATE & PORTABLE USER-FRIENDLY TOOLS TO MEASURE IOP VALUE COMPARABLE TO GOLDMANN APPLANATION TONOMETRY.

1. Lightweight (145gm w/o battery), portable and accuracy comparable with Goldmann Applanation Tonometry.
2. Can be used without anesthesia in OPD, OT, Patients on bed, ICU and Mass Glaucoma Screening Camp etc.
3. One year replacement warranty & conditional extended warranty against any manufacturing defect
4. Wireless data printing (Optional mini printer with built-in rechargeable battery).
5. Easy & user friendly operation, with L & R eye, IOP measurement in mmHg.




The SUOER Rebound Tonometer (model. SW-500) requires minimal operator training to assess accurate IOP, without the use of anesthesia, comparable well, with Goldman Applanation Tonometer. The Suoer Rebound Tonometer (Sw-500) is lightweight, portable and very user friendly for regular use in busy Eye clinics, IOP Screening in Eye Camps and also for multi clinic use.

Its unique self- calibration technology have made it as one of the Fast IOP measuring portable instrument with repeatability in accuracy and patient comfort in comparison to Goldmann Applanation Tonometry or Non-Contact Tonometry. Outstanding features such as; large patient memory with left and right eye datas and optional wireless mini pocket printer to take hard copy of measured data has super seeded suoer Rebound Tonometer (SW-500) over its all global rivals.

### PRINCIPLE

It follows a simple procedure where a very light disposable probe rebounds over the central cornea from a fixed distance (approx. 5mm) with practically invisible & painless impact, there by displaying the IOP in no time.

In this technology the motion parameters of the light probe is recorded through an induction based coil system. Then advanced algorithm combined with latest computer software, analyzes the deceleration of probe, its contact time and other parameters while touching the cornea. This deceleration and other rebound parameter of probe ultimately results the IOP. In brief, higher the IOP, faster the deceleration of probe and shorter the probe contact time.



The ratio of horizontal diameter of physiological cup to that of horizontal diameter of Optic Disc, is known to be Cup-Disc Ratio and should be less than 0.5. If it exceeds or there is a difference in the ratio between two eyes or if there is a progressive enlargement

, then GLAUCOMA may be suspected. One of the prime factors of GLAUCOMA is abnormal IOP. A perfect, fast and accurate tonometer that can read IOP values comparable to Goldmann Applanation tonometry is the **SUOER RE-BOUNDED TONOMETER : SW-500** a technological advancement of 21<sup>st</sup> Century.

6. Large memory to store, browse & print IOP data for Both L & R eyes of patients.
7. Presented with attractive Carry case for secure transportation.
8. Least Power Consumption ( 2xAA Dry Cells)



Cupping of Optic Disc

Portable NCT

Desktop NCT



# Technical Details of Suoer Rebound Tonometer : SW-500




Fig-1). Ultra light disposable probe used to measure IOP (resting on nozzle tip of SW-500) guarantees the safety from cross-infections.



Fig.-2). SW-500, probe-tip is held about 4-5 mm away from central cornea to touch it momentarily to read IOP, while Right or Left button is pushed.



Fig-3). Clear digital display shows 3 IOP readings with Avg. (P). R & L buttons for selecting Eye, Enter to memorize and print data from memory.




Fig.4). Portable, Handy and Shock-Proof packing in attractive carry case, for secure transportation.



Parts of SUOER : SW-500



Perfect Way of Using SW-500



Wireless Fast Thermal Printer




Fig-8.). Graph shows a direct comparison of readings between GAT & SW-500 RBT (studied on 75 patients). The standard deviation is found to be within  $\pm 2.00$ mmHg, a well accepted limit.




Fig-9). Graph shows the comfortability feelings of patients (75 patients tested with both GAT & SW-500 RBT). COMFORTABILITY PARAMETER:-  
 1.- Most Comfortable  
 2.- Comfortable  
 3.- Least Comfortable ;  
 4.- Irritating & uncomfortable .

Note: GAT: Goldmann Applanation Tonometry and RBT: Re-Bound Tonometry

TECHNICAL SPECIFICATION	
Measuring Range	: 3mmHg~70mmHg.
Measurement Accuracy	: $\pm 1.5$ mmHg (3mmHg $\leq$ Intra-Ocular pressures $\leq 25$ mmHg) : $\pm 2.5$ mmHg (25mmHg < Intra-Ocular pressures $\leq 70$ mmHg)
Weight	: 145 gm ( w/o Battery ) & 195 gm ( with batteries)
Power source	: 2 x AA batteries
Packing Dimension / Weight	: (L) 410 x (W) 320 x (H) 90mm / 2.50 Kg ( with Printer & Charger incl. Carry case)
Approval	: CE & ISO Certified.
Standard Package	: SW-500, box of probes (100 Nos.), One spare nozzle, Carry case with Key.
Optional Accessories	: Mini Wire less Fast thermal Printer with Built in Re-chargeable battery and one Paper Roll, Battery Charger, additional Probes ( box of 100 Nos.).

Exclusively Distributed by :



**Tianjin Suowei Electronic Technology Co. Ltd**  
 E-1-2F, No.6 Zhuyuan Road Huayuan Industrial Zone Nankai District  
 Tianjin-300384 P.R. China



**Mediscience Devices**

52/2B/2/1, Sarat Ghosh Garden Road, Kolkata - 700 031  
 Phone : (033) 2415 4090, Tele Fax : (033) 2415 1446  
 E-mail : mediscience@dataone.in / mediscience.calcutta@gmail.com  
[www.mediscience.tradeindia.com](http://www.mediscience.tradeindia.com) / [www.mediscience.in](http://www.mediscience.in)