

## Prima Mu OPH Specifications

45° inclinable binoculars. Optional 0-210° tiltable

3 step magnichanger 0.6x, 1.0x, 1.6x

Microscope carriers Single axis (Z) foot control

Eyepieces WF 10X/16 mm with eye guards, diopter adjustment  $\pm 5$  mm and diopter lock

Wide field focusable eyepieces 10x, 16mm FOV

Objective f-200mm. Optional f-175mm

Filters Green or Yellow

27W LED illumination

Optional Beam Splitters, HD Camera, Assistoscope

**Superb Red Reflex**



**Contact: 90380 81446**