



schiotz tonometer


schiotz A

schiotz B


7.5g



10g



5.5g


schiotz C

A masterpiece of precision mechanics! Highly precise eye tonometer according to Prof. Schiötz for measuring intra-ocular pressure. Its **reliability and durability** makes this quality measuring instrument indispensable in every ophthalmologist's practice and eye clinic.

- ▶ High quality agate bearing for an extremely long service life
- ▶ Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1 subdivisions
- ▶ Perfect reading of the scale with red pointer
- ▶ All vital parts are made of stainless steel, other parts are chrome-plated
- ▶ Complete in a black deluxe case with velvet-look inserts
- ▶ Supplied with three weights (5.5 g, 7.5 g, 10 g) and a conversion table

No. 5110	schiotz A Tonometer, straight scale, silver
No. 5111	schiotz A Tonometer, inclined scale, silver
No. 5112	schiotz A Tonometer, inclined scale, silver
No. 5113	schiotz A Tonometer, inclined scale, gold
No. 5114	schiotz C Tonometer, inclined scale, silver, with 4 weights

Please visit **riester.de** for more information