





FUJITRON

PORTABLE WIRELESS HANDY SLIT LAMP

PFSL-2X Pro. (with German Optics)





◆ EXCELLENT OPTICAL QUALITY

Especially fit for children, ward round, screening

◆ BIG SLIT ANGLE

Illumination angle max. up to be 60°, can get better slit view for cornea etc.

◆ WIDE VISUAL FIELD

Superb German Optics with 50% higher field of view than traditional design and a wider observation range

◆ ONE-HAND OPERATION

Easily adjust observation magnification during the inspection, help doctors quickly find the lesion in a larger field of vision and transfer to a higher magnification for better observation

◆ WARM LED

Warm light for better clinical observation of eye tissues and LED light source for longer and safer use

◆ SPECIFICATION

Magnification	10X & 16X
Working Distance	80mm(16X) 100mm(10X)
Field of View	ø10mm(16X) ø15mm(10X)
Interpupillary	49-75mm
Diopter Adjustment	±7D
Slit Size	0.1mm*12mm 0.2mm*8mm 0.8mm*12mm
Spot Size	ø1 / ø5 / ø12mm
Slit Angel	Horizontal, ±60°
Filter	Neutral Density Red-free Cobalt Blue
Light Resource	Warm LED
Illumination control	Continuously Adjustable
Battery Type	2 AA batteries
Operating Time	≥4 Hours
Dimensions	320*310*205mm
Net Weight	Approx. 835g



Manufactured for Fujitron Ophthalmics.
(An OEM brand of Mediscience Devices)

Authorised Distributors :



Mediscience Devices
(A Unit of Mediscience Services)

52/2B/2/1, Sarat Ghosh Garden Road (Dhakuria), Kolkata - 700 031
Mob. : +91 93306 76823 / 90380 81446
E-mail : mediscience.calcutta@gmail.com
Website : www.indiamart.com/mediscience-services