



SUOER SW-100 Digital Handheld Keratometer

Features:

Ultra small corneal curvature tester, is mainly used to measure the corneal curvature radius and diopter, wireless output print data.



Specifications:

Model	SW-100
Measuring Range	6.5mm~9.5mm: ±0.05mm: 0.01mm
Precision	± 0.05mm
Resolution of Curvature Radius of Cornea	0.01mm
Single Measuring Time	0.03s
Output	Wireless Infrared Thermal Printer
Can observe the eye directly through the screen.	
Weight	0.5Kg(with batteries)
Dimension	240mmx90mmx60mm
Power	500mW+15%