



MEDIVISION SOLID STATE BI-POLAR WETFIELD CAUTERY MBC-2020



Medivision Solid State Wet Field Bipolar Cautery

Salient Features

High Frequency Solid state circuitry with Pulse Wave Modulation (PWM) Technology.
Compact, light weight, elegant look and easy to use.
Reusable forceps, eraser (Optional) and cables.
Foot switch control.
Audio Visual indication (7-Segment Display).
Dial for variable power output (Continuous/Foot Control).

Technical Specifications

Supply voltage	: 220 - 230V AC
Frequency	: 50/60 Hz
Output power	: 10Watts (PWM Technology)
Power adjustment control	: Digitally controlled Knob
Display & Range	: 7-segment digital display / 0 – 10W
Dimension	: 200mm x 200mm x 70mm
Casing	: M.S powder coated
Weight	: 1 Kg approx.
Activation of HF power	: Foot switch (detachable)/Continuous Mode without Foot switch

Standard Accessories: Silicon Bipolar Cable (Autoclavable) 1set, Bipolarforceps 1 No., Foot Switch;1set