



MEDIVISION EXECUTIVE TRIAL LENS SET

A time tested quality product of Mediscience Devices, manufactured under intensive quality control in several stages. Each Lens is made of imported glass fitted in Colored Aluminum rims. A total of 190 Pcs of Lenses and 8 Pcs of accessories properly arranged in a slant tray fitted in a wooden box for easy handling. Medivision illuminated trial lens sets have established its brand image past 3 decades with 100s of optometrist and eye surgeons for its unparalleled quality and accuracy.

TECHNICAL SPECIFICATION

SPHERICAL CONCAVE (-Ve POWER/MYOPIC)			SPHERICAL CONVEX (+Ve POWER/HYPEROPIC)		
Diopter Range (-Ve)	Step	Total Lens	Diopter Range (+Ve)	Step	Total Lens
0.12	0.12	2	0.12	0.12	2
0.25 ~ 4D	0.25	32	0.25 ~ 4D	0.25	32
4.0 ~ 6D	0.50	8	4.0 ~ 6D	0.50	8
6.0 ~ 14D	1.0D	16	6.0 ~ 14D	1.0D	16
14.0 ~ 20D	2.0D	6	14.0 ~ 20D	2.0D	6
CYLINDER			CYLINDER		
0.12	0.12	2	0.12	0.12	2
0.25 ~ 3.50	0.25	14	0.25 ~ 3.50	0.25	14
3.50 ~ 6.0	0.50	10	3.50 ~ 6.0	0.50	10
ACCESSORIES	MADOX-1, SLIT-1, PIN HOLE-1, RED-1, GREEN-1, BLANK-2, ADULT TRIAL FRAME-1, LENS CLEANING CLOTH-1 Pc				
PRISMS	½, 1, 2, 3, 4, 5, 6, 8, 10 & 12 = 10 Pcs		TOTAL LENS - 190Pcs.		



Medivision Executive Trial lens Set.



Trial-Frame



Illuminated Dori Drum



Medivision Remote Control Mot. Vision Drum



Medivision LED Smart Chart



Ind. Medivision Streak Retinoscope



Keeler LED Ophthalmoscope & Retinoscope



Medivision Surgical Loop & Light

MEDIVISION DORI DRUM-- An economical Portable & robust illuminated reverse vision Drum(10 Ft) incorporated with 4 commonly used Language chart with C/E chart for illiterate. The respective chart is selected by a Pulling rope & roller arrangement.

MEDIVISION MOTORISED REMOTE CONTROLLED VISION DRUM-- An illuminated motorized reverse vision drum operated electronically through remote to move or select the required Snellen's chart.

MEDIVISION STREAK RETNOSCOPE: MVR-200-- An economical indigenous steak retinoscope available in "1.5V x 2" Dry Cell (Rs 2500/- + GST), as well as "3.6V" built-in rechargeable handle with charger (Rs 3700/- + GST)

MEDIVISION SURGICAL LOUPES & LIGHTS — Economical magnifying surgical binocular Loupes Available in 2.5x, 3.0x, 3.5x & 4.0x fitted over spectacle frame or head band. Also Available Wireless LED headlight (both diffused or spot light) which can be supplied in head band version or CLAMP FIT over Loupes with Spectacle Frame.


Mediscience Services
 52/2B/2/1, Sarat Ghosh Garden Road, Kolkata - 700 031 (W.B.)
 Phone : 9330676823 / 9038081446
 E-mail : mediscience.calcutta@gmail.com Web : www.mediscience.co.in