



**BIRO - TONOMETER - a new generation instrument from Germany - in quality that lasts !**

1.) Anodized scale mount

- highly resistant to sterilizing wear.

2.) Pointer and lever

- made of hardened stainless steel (Rockwell HRC 54), particularly appropriate for the manufacture of instruments; consequently no risk of pointer bending. Very attractive, gold plated design.

3.) Finger rest and transversal support


- black (chromium-plated), modern, pleasant design.

4.) 10 gram and 7.5 gram weight units

- with V-shaped keyway slit; for easy introduction of weight unit into tonometer.

5.) Metal case

- modern, finely padded steel case. For safe storage and protection of your BIRO-TONOMETER.



BIRO-TONOMETER, a new generation instrument from Germany - in quality that lasts! Manufactured in adherence to Professor Schioetz's theories. Suitable for measuring intraocular pressure. Complete with set of 3 weight units (5.5, 7.5 and 10 grams). Available with either straight scale layout or scale curved to an angle of 45°.

**Gerhard Biro**

Med.Apparate, Engeleswies 1  
 72393 Burladingen  
 E-mail : Gerhard.Biro@T-Online.de