

WIRELESS HANDHELD AUTOREFRACTOMETER INSTAREF : R20



Remidio Instaref R20

A handheld portable auto refractometer
for screening refractive errors



Accuracy and Mobility

Bringing in cutting-edge technology that travels with you.

With **Shack Hartmann aberrometry technology**, accurate reliable refractive values can be provided which saves you examination time.

Screen your patients wherever and however you see them

Examine patients anywhere without causing discomfort including camps, clinics and at home.

Fast and easy reading

In less than three seconds, you will have three quick and accurate readings and a final average value.





Fully Connected for the modern medical practice.

Reports can easily be printed using Bluetooth enabled wireless printers instantly.



Supported by clinical evidence

A clinical study of 122 patients against a predicate desktop AR device Topcon 8900 showed an acceptable equivalence within +/- 0.25D

Docking station for easy charging





Product Specifications

Model / Part no.	Instaref R20 / 1000010021
Methodology	Shack Hartmann wavefront
Sphere range	-10.0 D to +10.0D
Cylinder range	-5.0 D to +5.0D
Axis of astigmatism range	1-180°
Connection	USB, Bluetooth
Use mode	Handheld
Weight	
• Device	450g
• Device with docking station	750g



Marketed by :

MEDISCIENCE

52/2B/2/1, SARAT GHOSH ROAD, KOLKATA - 700031(WB), INDIA.

E-mail : mediscience.calcutta@gmail.com; Phone : 9038081446 / 9354127547 / 9330676823

