



Mediscience Devices
(An Unit of Mediscience Services)

BALIWALA & HOMI NON-ILLUMINATED TRIAL LENS SET

The Baliwala & Homi Non-Illuminated Trial Lens Set is a high-quality ophthalmic diagnostic kit designed for precise and reliable subjective refraction. The set includes a comprehensive range of premium optical glass trial lenses housed in a durable non-illuminated carrying case for organized storage and easy access. Supplied with a precision adjustable trial frame, this trial lens set is ideal for ophthalmologists, optometrists, eye hospitals, optical stores, vision care centers, and medical colleges.

KEY FEATURES

- Premium Non-Illuminated Trial Lens Set
- High-Quality Optical Glass Trial Lenses
- Anti-Reflective (AR) Coated Lenses
- Precision Adjustable Trial Frame Included
- Complete Sphere, Cylinder & Prism Lens Set
- Plus & Minus Sphere Lenses
- Plus & Minus Cylinder Lenses
- Maddox Rod, Pinhole Disc & Occluder Included
- Cross Cylinder & Frosted Lens Included
- Color-Coded Lens Identification
- Durable Aluminum Carrying Case
- Easy Lens Organization & Quick Access
- Lightweight & Portable Design
- Ideal for Ophthalmologists, Optometrists & Eye Clinics



SPECIFICATIONS

Parameter	Specification
Product Type	Non-Illuminated Trial Lens Set
Brand	Baliwala & Homi
Lens Material	High-Quality Optical Glass
Lens Coating	Anti-Reflective (AR) Coated
Sphere Lenses	Plus & Minus Set
Cylinder Lenses	Plus & Minus Set
Prism Lenses	Included
Trial Frame	Precision Adjustable Metal Trial Frame
Accessories	Maddox Rod, Pinhole Disc, Occluder, Frosted Lens & Cross Cylinder
Lens Identification	Color Coded
Storage Case	Durable Aluminum Carrying Case
Application	Subjective Refraction & Vision Testing
Suitable For	Optical Shops, Eye Hospitals & Ophthalmology Clinics

DISTRIBUTED BY:
MEDISCIENCE DEVICES

9330676823
mediscience.calcutta@gmail.com
www.mediscience.co.in
2/1, 52/2 B, Sarat Ghosh Garden Rd, Dhakuria, Kasba,
Kolkata, West Bengal 700031